

Chirurgische Nadeln

Surgical Needles


CEATEC® MEDIZINTECHNIK

Nadelformen und Zeichenerklärung:
Needle Shapes and sign explanation:

Nadelform:	***Needle Shape:***
3/8 kreisförmig	*3/8 circle*
1/2 kreisförmig	*1/2 circle*
halb-gebogen	half curved
gerade	*straight*

Nadelart:		***Needle Type:***
dreikantig, innen-liegende Schneide		three-sided, *cutting edge inside*
runde Spitze		*round/tapered point*

Nadelöhr:		***Needle Eye:***
gewöhnliches Öhr		regular eye
Federöhr		*spiral eye*




= gestreckte Länge
= *straightened length*

Verpackungsseinheit (VE): Päckchen á 12 Stück im Minigrip-Beutel
Packing Unit (PU): Package of 12 pieces in MInigrip-Bag

B


gebogen 3/8 Kreis, VE = 12 Stück
curved 3/8 circle, SU = 12 pieces



				Figur *Figure*	gestreckte Länge *Straightened Length*
AE 112G01	AE 113G01	AE 114G01	AE 115G01	1	1,40 x 80 mm
AE 112G.02	AE 113G02	AE 114G02	AE 115G02	2	1,40 x 73 mm
AE 112G03	AE 113G03	AE 114G03	AE 115G03	3	1,30 x 68 mm
AE 112G04	AE 113G04	AE 114G04	AE 115G04	4	1,30 x 63 mm
AE 112G05	AE 113G05	AE 114G05	AE 115G05	5	1,20 x 58 mm
AE 112G06	AE 113G06	AE 114G06	AE 115G06	6	1,10 x 54 mm
AE 112G07	AE 113G07	AE 114G07	AE 115G07	7	1,10 x 50 mm
AE 112G08	AE 113G08	AE 114G08	AE 115G08	8	1,10 x 46 mm
AE 112G09	AE 113G09	AE 114G09	AE 115G09	9	1,00 x 42 mm
AE 112G10	AE 113G10	AE 114G10	AE 115G10	10	0,90 x 38 mm
AE 112G11	AE 113G11	AE 114G11	AE 115G11	11	0,80 x 34 mm
AE 112G12	AE 113G12	AE 114G12	AE 115G12	12	0,70 x 30 mm
AE 112G13	AE 113G13	AE 114G13	AE 115G13	13	0,70 x 27 mm
AE 112G14	AE 113G14	AE 114G14	AE 115G14	14	0,70 x 23 mm
AE 112G15	AE 113G15	AE 114G15	AE 115G15	15	0,60 x 20 mm
AE 112G16	AE 113G16	AE 114G16	AE 115G16	16	0,60 x 18 mm

G


gebogen 1/2 Kreis, VE = 12 Stück
curved 1/2 circle, SU = 12 pieces



				Figur *Figure*	gestreckte Länge *Straightened Length*
AE 122G01	AE 123G01	AE 124G01	AE 125G01	1	1,40 x 80 mm
AE 122G.02	AE 123G02	AE 124G02	AE 125G02	2	1,40 x 73 mm
AE 122G03	AE 123G03	AE 124G03	AE 125G03	3	1,30 x 68 mm
AE 122G04	AE 123G04	AE 124G04	AE 125G04	4	1,30 x 63 mm
AE 122G05	AE 123G05	AE 124G05	AE 125G05	5	1,20 x 58 mm
AE 122G06	AE 123G06	AE 124G06	AE 125G06	6	1,10 x 54 mm
AE 122G07	AE 123G07	AE 124G07	AE 125G07	7	1,10 x 50 mm
AE 122G08	AE 123G08	AE 124G08	AE 125G08	8	1,10 x 46 mm
AE 122G09	AE 123G09	AE 124G09	AE 125G09	9	1,00 x 42 mm
AE 122G10	AE 123G10	AE 124G10	AE 125G10	10	0,90 x 38 mm
AE 122G11	AE 123G11	AE 124G11	AE 125G11	11	0,80 x 34 mm
AE 122G12	AE 123G12	AE 124G12	AE 125G12	12	0,70 x 30 mm
AE 122G13	AE 123G13	AE 114G13	AE 125G13	13	0,70 x 27 mm
AE 122G14	AE 123G14	AE 124G14	AE 125G14	14	0,70 x 23 mm
AE 122G15	AE 123G15	AE 124G15	AE 125G15	15	0,60 x 20 mm
AE 122G16	AE 123G16	AE 124G16	AE 125G16	16	0,60 x 18 mm

Ga


gebogen 1/2 Kreis, VE = 12 Stück
curved 1/2 circle, SU = 12 pieces



				Figur *Figure*	gestreckte Länge *Straightened Length*
AE 132G01	AE 133G01	AE 134G01	AE 135G01	1	1,60 x 65 mm
AE 132G.02	AE 133G02	AE 134G02	AE 135G02	2	1,60 x 60 mm
AE 132G03	AE 133G03	AE 134G03	AE 135G03	3	1,40 x 52 mm
AE 132G04	AE 133G04	AE 134G04	AE 135G04	4	1,40 x 45 mm
AE 132G05	AE 133G05	AE 134G05	AE 135G05	5	1,30 x 37 mm
AE 132G06	AE 133G06	AE 134G06	AE 135G06	6	1,10 x 32 mm
AE 132G07	AE 133G07	AE 134G07	AE 135G07	7	1,00 x 26 mm
AE 132G08	AE 133G08	AE 134G08	AE 135G08	8	1,00 x 21 mm
AE 132G09	AE 133G09	AE 134G09	AE 135G09	9	0,90 x 17 mm

E


gebogen 1/2 Kreis, VE = 12 Stück
curved 1/2 circle, SU = 12 pieces





		Figur *Figure*	gestreckte Länge *Straightened Length*
AE 154G30	AE 155G30	000	1,10 x 72 mm
AE 154G20	AE 155G20	00	1,10 x 63 mm
AE 154G00	AE 155G00	0	1,00 x 54 mm
AE 154G01	AE 155G01	1	1,00 x 45 mm
AE 154G02	AE 155G02	2	0,90 x 36 mm
AE 154G03	AE 155G03	3	0,80 x 29 mm
AE 154G04	AE 155G04	4	0,70 x 23 mm
AE 154G05	AE 155G05	5	0,60 x 20 mm

Keith

gerade, VE = 12 Stück
straight, SU = 12 pieces





				**Figur** ***Figure***	**gestreckte Länge** ***Straightened Length***
AE 162G01	AE 163G01	AE 164G01	AE 165G01	1	0,80 x 75 mm
AE 162G.02	AE 163G02	AE 164G02	AE 165G02	2	0,80 x 70 mm
AE 162G03	AE 163G03	AE 164G03	AE 165G03	3	0,80 x 65 mm
AE 162G04	AE 163G04	AE 164G04	AE 165G04	4	0,70 x 60 mm
AE 162G05	AE 163G05	AE 164G05	AE 165G05	5	0,70 x 55 mm
AE 162G06	AE 163G06	AE 164G06	AE 165G06	6	0,70 x 50 mm

Darmnadel
Intestinal Needle

Pa


halb-gebogen, VE = 12 Stück
half-curved, SU = 12 pieces



		**Figur** ***Figure***	**gestreckte Länge** ***Straightened Length***
AE 170G20	AE 171G20	00	0,70 x 68 mm
AE 170G00	AE 171G00	0	0,70 x 60 mm
AE 170G01	AE 171G01	1	0,70 x 52 mm
AE 170G02	AE 171G02	2	0,70 x 42 mm
AE 170G03	AE 171G03	3	0,60 x 36 mm
AE 170G04	AE 171G04	4	0,60 x 30 mm
AE 170G05	AE 171G05	5	0,60 x 24 mm

Darmnadel
Intestinal Needle


Pb gebogen 1/2 Kreis, VE = 12 Stück
curved 1/2 circle, SU = 12 pieces



		Figur *Figure*	gestreckte Länge *Straightened Length*
AE 174G20	AE 175G20	00	0,70 x 68 mm
AE 174G00	AE 175G00	0	0,70 x 60 mm
AE 174G01	AE 175G01	1	0,70 x 52 mm
AE 174G02	AE 175G02	2	0,70 x 42 mm
AE 174G03	AE 175G03	3	0,60 x 36 mm
AE 174G04	AE 175G04	4	0,60 x 30 mm
AE 174G05	AE 175G05	5	0,60 x 24 mm

Darmnadel
Intestinal Needle


Pc gerade, VE = 12 Stück
straight, SU = 12 pieces



		Figur *Figure*	gestreckte Länge *Straightened Length*
AE 180G20	AE 181G20	00	0,70 x 68 mm
AE 180G00	AE 181G00	0	0,70 x 60 mm
AE 180G01	AE 181G01	1	0,70 x 52 mm
AE 180G02	AE 181G02	2	0,70 x 42 mm
AE 180G03	AE 181G03	3	0,60 x 36 mm
AE 180G04	AE 181G04	4	0,60 x 30 mm
AE 180G05	AE 181G05	5	0,60 x 24 mm

Darmnadel
Intestinal Needle


Pd gebogen 3/8 Kreis, VE = 12 Stück
curved 3/8 circle, SU = 12 pieces



		Figur *Figure*	gestreckte Länge *Straightened Length*
AE 184G20	AE 185G20	00	0,70 x 68 mm
AE 184G00	AE 185G00	0	0,70 x 60 mm
AE 184G01	AE 185G01	1	0,70 x 52 mm
AE 184G02	AE 185G02	2	0,70 x 42 mm
AE 184G03	AE 185G03	3	0,60 x 36 mm
AE 184G04	AE 185G04	4	0,60 x 30 mm
AE 184G05	AE 185G05	5	0,60 x 24 mm

Augennadel
Eye Needle


Oa halb-gebogen, VE = 12 Stück
half-curved, SU = 12 pieces



		Figur *Figure*	gestreckte Länge *Straightened Length*
AE 192G00	AE 193G00	0	0,60 x 32 mm
AE 192G01	AE 193G01	1	0,60 x 28 mm
AE 192G02	AE 193G02	2	0,60 x 24 mm
AE 192G03	AE 193G03	3	0,60 x 20 mm
AE 192G04	AE 193G04	4	0,60 x 15 mm
AE 192G05	AE 193G05	5	0,60 x 12 mm

Augennadel
Eye Needle

Ob gebogen 3/8 Kreis, VE = 12 Stück
curved 3/8 circle, SU = 12 pieces




		Figur *Figure*	gestreckte Länge *Straightened Length*
AE 196G00	AE 197G00	0	0,60 x 32 mm
AE 196G01	AE 197G01	1	0,60 x 28 mm
AE 196G02	AE 197G02	2	0,60 x 24 mm
AE 196G03	AE 197G03	3	0,60 x 20 mm
AE 196G04	AE 197G04	4	0,60 x 15 mm
AE 196G05	AE 197G05	5	0,60 x 12 mm

Augennadel
Eye Needle

Og

gebogen 1/2 Kreis, VE = 12 Stück
curved 1/2 circle, SU = 12 pieces




		**Figur** ***Figure***	**gestreckte Länge** ***Straightened Length***
AE 200G01	AE 201G01	1	0,60 x 28 mm
AE 200G02	AE 201G02	2	0,60 x 24 mm
AE 200G03	AE 201G03	3	0,60 x 20 mm
AE 200G04	AE 201G04	4	0,60 x 15 mm
AE 200G05	AE 201G05	5	0,60 x 12 mm

Gefäßnadel
Vascular Needle

Kalt

gebogen 3/8 Kreis, VE = 12 Stück
curved 3/8 circle, SU = 12 pieces




	**Figur** ***Figure***
AE 204G20	00
AE 204G00	0
AE 204G01	1
AE 204G02	2
AE 204G03	3

Gefäßnadel
Vascular Needle

Payr

gerade, VE = 12 Stück
straight, SU = 12 pieces





	**Figur** ***Figure***
AE 205G01	1
AE 205G02	2
AE 205G03	3
AE 205G04	4
AE 205G05	5

Gefäßnadel
Vascular Needle

Payr gebogen 3/8 Kreis, VE = 12 Stück
curved 3/8 circle, SU = 12 pieces



	**Figur** ***Figure***
AE 206G20	00
AE 206G00	0
AE 206G01	1
AE 206G02	2
AE 206G03	3






Abmessung in mm
Sizes in mm
l x b x h


SR 300.50	55 x 25 x 4,5
SR 300.65	65 x 42 x 8
SR 300.75	75 x 25 x 4,5



Abmessung in mm
Sizes in mm
l x b x h


SR 302.50	55 x 25 x 4,5
SR 302.65	65 x 42 x 8
SR 302.75	75 x 25 x 4,5





Abmessung in mm
Sizes in mm
Ø x h

SR 309.65	65 x 15



Innendose allein,
Inner case only

Abmessung in mm
Sizes in mm
Ø x h

SR 310.65	62 x 13



SAUERBRUCH

Abmessung in mm
Sizes in mm
l x b x h


SR 162.20	210 x 125 x 17



SAUERBRUCH

Abmessung in mm
Sizes in mm
l x b x h

SR 164.20	210 x 125 x 17





Abmessung in mm
Sizes in mm
l x b x h

SR 335.01 150 x 95 x 13



Abmessung in mm
Sizes in mm
l x b x h

SR 335.02 150 x 95 x 13



Abmessung in mm
Sizes in mm
l x b x h

SR 335.03 150 x 95 x 13



Abmessung in mm
Sizes in mm
l x b x h

SR 335.05 150 x 95 x 13




Abmessung in mm
Sizes in mm
l x b x h

SR 337.01 150 x 95 x 13



Abmessung in mm
Sizes in mm
l x b x h

SR 337.02 150 x 95 x 13



Abmessung in mm
Sizes in mm
l x b x h

SR 340.02 150 x 95 x 11

CEATEC Medizintechnik GmbH
Untere Hauptstrasse 6
D - 78573 Wurmlingen

web: www.ceatec.de
e-mail: info@ceatec.de
Tel.: +49-7461-9348-0
Fax: +49-7461-9348-19

05/2013

Technische Änderungen vorbehalten.
Technical modifications subject to change without notice.

Con riserve di modifiche tecniche.
Reservado el derecho de efectuar modificaciones técnicas.